

Anti-Crip1 Polyclonal Antibody For Wb And Elisa - P50238

Item Number: CM0030453

Introduction

Rabbit polyclonal antibody targeting human CRIP1 (cysteine-rich protein 1), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Recombinant full-length immunogen. For zinc transport research.

[Learn More](#)

Parameter	Value
Product Name	CRIP1 Rabbit pAb
Remarks/Alias	CRHP; CRIP; CRP1; CRP-1; CRIP1
Species	Human
GeneID	1396
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-77 of human CRIP1 (NP_001302.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50238
Protein Weight	9 kDa
Shipping	Ice bag