

Anti-Crot Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9Ukg9

Item Number: CM0027774

Introduction

Anti-CROT Rabbit Polyclonal Antibody validated for WB, IP, ELISA in Human, Mouse, Rat.

Recognizes peroxisomal carnitine O-octanoyltransferase (70 kDa), a key enzyme in fatty acid beta-oxidation.

[Learn More](#)

Parameter	Value
Product Name	CROT Rabbit pAb
Remarks/Alias	COT; CROT
Species	Human
GeneID (Human)	54677
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-87 of human CROT (NP_001230674.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKG9
Protein Weight	70kDa
Shipping	Ice bag