


Anti-Cryab Rabbit Monoclonal Antibody For Ihc-P, If/Icc, Elisa - P02511

Item Number: CM0004790



Introduction

Rabbit monoclonal antibody against human CRYAB (Alpha-crystallin B chain), validated for IHC-P, IF/ICC, ELISA. Cross-reacts with rat. Ideal for studying chaperone activity and lens biology in research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	CRYAB Rabbit mAb
Remarks/Alias	MFM2, CRYA2, CTPP2, HSPB5, CMD1II, CTRCT16, HEL-S-101, CRYAB
Species	Human
GeneID (human)	1410
GeneID	1410
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 76-175 of human CRYAB (P02511).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Rat
SWISS	P02511
Protein Weight	20kDa
Shipping	Ice bag