

Anti-Cryz Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q08257

Item Number: CM0018014



Introduction

High-quality rabbit polyclonal antibody targeting human CRYZ (zeta-crystallin/NADPH:quinone reductase). Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studies on mRNA stability, quinone detoxification, and oxidative stress.

[Learn More](#)

Parameter	Specification
Product Name	CRYZ Rabbit pAb
Remarks/Alias	CRYZ
Species	Human
GeneID (Human)	1429
GeneID	1429
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-329 of human CRYZ (NP_001123514.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08257
Protein Weight	35kDa
Shipping	Ice bag