

Anti-Cs Polyclonal Antibody For Wb, If/Icc, Elisa - 075390

Item Number: CM0025551



Introduction

Rabbit polyclonal antibody against human CS (citrate synthase, mitochondrial). Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for studying the Krebs cycle and mitochondrial metabolism.

[Learn More](#)

Parameter	Value
Product Name	CS Rabbit pAb
Remarks/Alias	CS
Species	Human
Gene ID (Human)	1431
Gene ID	1431
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 167-466 of human CS (NP_004068.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75390
Protein Weight	52kDa
Shipping	Ice bag