

Anti-Csad Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y600

Item Number: CM0017207



Introduction

Rabbit polyclonal antibody against human CSAD (cysteine sulfinic acid decarboxylase), for WB and ELISA. Cross-reacts with mouse and rat. Ideal for taurine biosynthesis research and liver/brain studies.

[Learn More](#)

Parameter	Value
Product Name	CSAD Rabbit pAb
Remarks/Alias	CSD; PCAP; CSAD
Species	Human
GeneID (human)	51380
GeneID	51380
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 371-520 of human CSAD (NP_057073.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y600
Protein Weight	55kDa
Shipping	Ice bag