

Anti-Csde1 Rabbit Monoclonal Antibody For Wb, Elisa - 075534

Item Number: CM0002503



Introduction

High-quality rabbit monoclonal antibody against CSDE1 (Cold shock domain-containing protein E1), validated for WB and ELISA. Ideal for studies on mRNA decay and translational regulation, with cross-reactivity to human and mouse.

[Learn More](#)

Parameter	Value
Product Name	CSDE1 Rabbit mAb
Remarks/Alias	UNR; D1S155E; CSDE1
Species	Human
GeneID (Human)	7812
GeneID	7812
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CSDE1 (NP_009089.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	075534
Protein Weight	89kDa
Shipping	Ice bag