


Anti-Csnk1G2 Rabbit Polyclonal Antibody For Wb, Elisa - P78368

Item Number: CM0025091



Introduction

High-quality rabbit polyclonal antibody against Casein kinase I isoform gamma-2 (CSNK1G2), validated for WB and ELISA. Detects human, mouse, and rat CSNK1G2, ideal for Wnt signaling and kinase research studies.

[Learn More](#)

Parameter	Value
Product Name	CSNK1G2 Rabbit pAb
Remarks/Alias	CK1g2; CSNK1G2
Species	Human
Gene ID (Human)	1455
Gene ID	1455
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 316-415 of human CSNK1G2 (NP_001310.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78368
Protein Weight	47kDa
Shipping	Ice bag