


Anti-Cstf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q05048

Item Number: CM0024955



Introduction

High-quality rabbit polyclonal antibody targeting human CSTF1, suitable for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying pre-mRNA 3'-end processing and polyadenylation.

[Learn More](#)

Parameter	Value
Product Name	CSTF1 Rabbit pAb
Remarks/Alias	CstF-50; CstFp50; CSTF1
Species	Human
GeneID (Human)	1477
GeneID	1477
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human CSTF1 (NP_001315.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05048
Protein Weight	48kDa
Shipping	Ice bag