

Anti-Ctb4A Polyclonal Antibody For Wb, Elisa - Q0Jf76

Item Number: CM0022646

Introduction

Rabbit polyclonal antibody targeting CTB4a (LRR-PK1) from Oryza sativa, suitable for WB and ELISA. Supports research on plant receptor-like kinases and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	CTB4a Rabbit pAb
Remarks/Alias	LRR-PK1; CTB4a
Species	Oryza sativa
GeneID	4335004
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-152 of Oryza sativa CTB4a. (Q0JF76)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q0JF76
Protein Weight	120kDa
Shipping	Ice bag