

Anti-Ctsf Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ubx1

Item Number: CM0016727



Introduction

Rabbit polyclonal antibody against human CTSF (cathepsin F), suitable for Western blot, IHC-P, and ELISA. Designed for research on lysosomal thiol protease roles in protein degradation and tumor invasion.

[Learn More](#)

Parameter	Value
Product Name	CTS F Rabbit pAb
Remarks/Alias	CATSF; CLN13; CTSF
Species	Human
Gene ID (Human)	8722
Gene ID	8722
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-484 of human CTSF (NP_003784.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9UBX1
Protein Weight	53kDa
Shipping	Ice bag