

Anti-Cullin 3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13618

Item Number: CM0020254



Introduction

Rabbit polyclonal antibody against human Cullin 3 (CUL-3, PHA2E), validated for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. For ubiquitination and proteasome pathway research.

[Learn More](#)

Parameter	Value
Product Name	Cullin 3 Rabbit pAb
Remarks/Alias	CUL-3; PHA2E; NEDAUS; Cullin 3
Species	Human
Gene ID (Human)	8452
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 569-768 of human Cullin 3 (NP_003581.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13618
Protein Weight	89kDa
Shipping	Ice bag