

Anti-Cullin-4A Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13619

Item Number: CM0009111



Introduction

Rabbit monoclonal antibody targeting human Cullin-4A (CUL4A), validated for WB, IHC-P, and ELISA. Recognizes an 88 kDa protein; suitable for ubiquitination and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	Cullin 4A Rabbit mAb
Remarks/Alias	CUL4A; cullin-4A; Cullin 4A
Species	Human
GeneID (Human)	8451
GeneID	8451
Immunogen	Synthetic peptide corresponding to amino acids 660-759 of human Cullin 4A (Q13619)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q13619
Protein Weight	88kDa
Shipping	Ice bag