

Anti-Cxcl9/Mig Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q07325

Item Number: CM0010867

Introduction

High-quality Anti-CXCL9/MIG Rabbit Monoclonal Antibody validated for WB, IF/ICC, ELISA. Specific to human CXCL9, a chemokine that directs T-cell migration via CXCR3. Ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	CXCL9/MIG Rabbit mAb
Remarks/Alias	CMK; MIG; Humig; SCYB9; crg-10; CXCL9
Species	Human
Gene ID (Human)	4283
Gene ID	4283
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-125 of human CXCL9/MIG(NP_002407.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q07325
Protein Weight	14kDa
Shipping	Ice bag