


Anti-Cxxc1 Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9P0U4

Item Number: CM0013603



Introduction

Rabbit polyclonal antibody targeting CXXC1 (CFP1/CpG-binding protein) for WB, IHC-P, IP, and ELISA. Validated in human and mouse samples; essential for transcription regulation and chromatin biology research.

[Learn More](#)

Parameter	Value
Product Name	CXXC1 Rabbit pAb
Remarks/Alias	CFP1; CGBP; SPP1; PCCX1; PHF18; hCGBP; ZCGPC1; HsT2645; 2410002I16Rik; 5830420C16Rik; CXXC1
Species	Human
GeneID (human)	30827
GeneID	30827
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human CXXC1 (NP_001095124.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P0U4
Protein Weight	76 kDa
Shipping	Ice bag