

Anti-Cyclin E1 Rabbit Polyclonal Antibody For Wb And Elisa - P24864

Item Number: CM0015962



Introduction

Rabbit polyclonal antibody targeting human Cyclin E1, validated for Western blot (WB) and ELISA. Shows cross-reactivity with mouse and rat, making it suitable for cell cycle and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Cyclin E1 Rabbit pAb
Remarks/Alias	CCNE; pCCNE1; Cyclin E1
Species	Human
GeneID (human)	898
GeneID	898
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-400 of human Cyclin E1 (NP_001229.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P24864
Protein Weight	47kDa
Shipping	Ice bag