

Anti-Cyclin K Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O75909

Item Number: CM0015740



Introduction

Cyclin K Rabbit Polyclonal Antibody (O75909) validated for WB, IHC-P, ELISA in human, mouse, rat. Regulatory subunit of CDKs; functions in transcription regulation via RNA Pol II CTD phosphorylation.

[Learn More](#)

Parameter	Value
Product Name	Cyclin K Rabbit pAb
Remarks/Alias	CPR4; IDDHDF; Cyclin K
Species	Human
Gene ID	8812 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human Cyclin K (NP_001092872.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75909
Protein Weight	64kDa
Shipping	Ice bag