

Anti-Cyclophilin 40 Polyclonal Antibody For Wb, Elisa - Q08752

Item Number: CM0021797

Introduction

Rabbit polyclonal antibody targeting human Cyclophilin 40 (PPID, CYP-40). Validated for WB and ELISA, cross-reactive with rat. Suitable for research on protein folding, HSP90 co-chaperone complexes, apoptosis, and viral replication.

[Learn More](#)

Parameter	Value
Product Name	Cyclophilin 40 Rabbit pAb
Remarks/Alias	CYPD; CYP-40; Cyclophilin 40
Species	Human
Gene ID (Human)	5481
Gene ID	5481
Immunogen	Recombinant protein corresponding to a sequence within amino acids 161-370 of human Cyclophilin 40 (NP_005029.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q08752
Protein Weight	41kDa
Shipping	Ice bag