

Anti-Cyp1A1 Polyclonal Antibody For Wb, If/Icc, Elisa - P04798

Item Number: CM0020094



Introduction

Anti-CYP1A1 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat CYP1A1. Recombinant immunogen (aa 200-460). Ideal for cytochrome P450 metabolism and toxicology research.

[Learn More](#)

Parameter	Value
Product Name	CYP1A1 Rabbit pAb
Remarks/Alias	AHH; AHRR; CP11; CYP1; CYP1A1; P1-450; P450-C; P450DX; CYP1A1
Species	Human
Gene ID (Human)	1543
Gene ID	1543
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-460 of human CYP1A1 (NP_000490.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04798
Protein Weight	58kDa
Shipping	Ice bag