


Anti-Cyp1A1 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P04798

Item Number: CM0023857



Introduction

Rabbit polyclonal antibody against human CYP1A1, validated for WB, IF-P, IHC-P, and ELISA. Useful for studying drug metabolism, hormone regulation, and lipid signaling. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	CYP1A1 Rabbit pAb
Remarks/Alias	AHH; AHRR; CP11; CYP1; CYP1A1; P1-450; P450-C; P450DX
Species	Human
GeneID (Human)	1543
GeneID	1543
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 261-360 of human CYP1A1 (NP_000490.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04798
Protein Weight	58kDa
Shipping	Ice bag