

Anti-Cyp26A1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - O43174

Item Number: CM0008852

Introduction

High-quality rabbit monoclonal antibody against human CYP26A1, for WB, IF-P, IHC-P, ELISA. Detects human, mouse, and rat CYP26A1. Target retinoic acid metabolism and signaling in various tissues.

[Learn More](#)

Parameter	Specification
Product Name	CYP26A1 Rabbit mAb
Remarks/Alias	CP26; CYP26; P450RAI; P450RAI1; CYP26A1
Species	Human
GeneID (Human)	1592
GeneID	1592
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 201-450 of human CYP26A1 (O43174).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43174
Protein Weight	56 kDa
Shipping	Ice bag