

Anti-Cyp39A1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Nyl5

Item Number: CM0012562

Introduction

Rabbit monoclonal antibody targeting CYP39A1, validated for WB, IHC-P, ELISA with cross-reactivity to Mouse. Supports research in cholesterol metabolism and IVD assay development. 54kDa predicted.

[Learn More](#)

Parameter	Value
Product Name	CYP39A1 Rabbit mAb
Species	Human
Gene ID (Human)	51302
Gene ID	51302
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 301-400 of human CYP39A1. (NP_057677.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	Q9NYL5
Protein Weight	54kDa
Shipping	Ice bag