

Anti-Cytochrome P450 2E1 (Cyp2E1) Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P05181

Item Number: CM0007237



Introduction

Rabbit monoclonal antibody against human CYP2E1, validated for WB, IF, IHC, and ELISA. Detects human, mouse, rat CYP2E1. Ideal for research in drug metabolism and toxicology.

[Learn More](#)

Parameter	Value
Product Name	Cytochrome P450 2E1 (CYP2E1) Rabbit mAb
Remarks/Alias	CPE1; CYP2E; P450-J; P450C2E; Cytochrome P450 2E1 (CYP2E1)
Species	Human
Gene ID	1571
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 234-493 of human Cytochrome P450 2E1 (CYP2E1) (NP_000764.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05181
Protein Weight	57kDa
Shipping	Ice bag