

Anti-Cytochrome P450 4A (Cyp4A11) Monoclonal Antibody For Wb, Ihc-P, Elisa - Q02928

Item Number: CM0002834

Introduction

Validated rabbit monoclonal antibody targeting human CYP4A11, suitable for WB, IHC-P, and ELISA. Enables research on fatty acid omega-hydroxylation, arachidonic acid metabolism, and 20-HETE signaling. Detects human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	Cytochrome P450 4A (CYP4A11) Rabbit mAb
Remarks/Alias	CP4Y; CYP4A2; CYP4AII; CYPIVA11; Cytochrome P450 4A (CYP4A11)
Species	Human
GeneID (Human)	1579
GeneID	1579
Immunogen	Synthetic peptide corresponding to amino acids 50-150 of human CYP4A11 (Q02928)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02928
Protein Weight	59kDa
Shipping	Ice bag