

## Anti-Cytokeratin 18 (Krt18) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P05783

Item Number: CM0017592

### Introduction

Rabbit polyclonal antibody against human Cytokeratin 18 (KRT18), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for epithelial cell biology and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                                                |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Cytokeratin 18 (KRT18) Rabbit pAb                                                                                                    |
| Remarks/Alias    | K18; CK-18; CYK18; Cytokeratin 18 (KRT18)                                                                                            |
| Species          | Human                                                                                                                                |
| Gene ID (Human)  | 3875                                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 127-430 of human Cytokeratin 18 (KRT18) (NP_000215.1). |
| Source           | Rabbit                                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                                |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                    |
| SWISS            | P05783                                                                                                                               |
| Protein Weight   | 48kDa                                                                                                                                |
| Shipping         | Ice bag                                                                                                                              |