


Anti-Cytokeratin 19 (Krt19) Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08727

Item Number: CM0016636



Introduction

Rabbit polyclonal antibody against human Cytokeratin 19 (KRT19/CK19), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for epithelial biology and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	Cytokeratin 19 (KRT19) Rabbit pAb
Remarks/Alias	K19; CK19; K1CS; Cytokeratin 19 (KRT19)
Species	Human
GeneID (human)	3880
GeneID	3880
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 210-340 of human Cytokeratin 19 (KRT19) (NP_002267.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08727
Protein Weight	44kDa
Shipping	Ice bag