

Anti-Cytokeratin 20 (Krt20) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P35900

Item Number: CM0013322

Introduction

Rabbit polyclonal antibody against Cytokeratin 20 (KRT20), suitable for WB, IF/ICC, ELISA. Reacts with human and mouse samples. Ideal for studying keratin filament organization, mucin secretion, and intestinal epithelium biology. Part of CamelBio's IVD raw material portfolio for diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 20 (KRT20) Rabbit pAb
Remarks/Alias	K20; CD20; CK20; CK-20; KRT21; Cytokeratin 20 (KRT20)
Species	Human
Gene ID	54474 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 245-424 of human Cytokeratin 20 (KRT20) (NP_061883.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P35900
Protein Weight	48kDa
Shipping	Ice bag