

Anti-Cytokeratin 2E (Krt2) Rabbit Polyclonal Antibody For Wb And Elisa - P35908

Item Number: CM0016174

Introduction

Rabbit polyclonal antibody targeting human Cytokeratin 2e (KRT2), validated for Western blot and ELISA applications. Cross-reacts with mouse and rat. Suitable for studies on epidermal differentiation, keratin intermediate filaments, and skin barrier function.

[Learn More](#)

Parameter	Value
Product Name	Cytokeratin 2e (KRT2) Rabbit pAb
Remarks/Alias	K2e; KRTE; CK-2e; KRT2A; KRT2E; Cytokeratin 2e (KRT2)
Species	Human
GeneID (Human)	3849
GeneID	3849
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human Cytokeratin 2e (KRT2) (NP_000414.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35908
Protein Weight	65kDa
Shipping	Ice bag