


## Anti-Cytokeratin 5 Monoclonal Antibody For Wb, Ihc-P, Elisa - P13647

Item Number: CM0002494



### Introduction

Rabbit monoclonal antibody against human Cytokeratin 5 (KRT5), suitable for Western blotting, immunohistochemistry (paraffin), and ELISA. Reacts with human, mouse, and rat samples. Ideal for epidermal and keratinocyte research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | Cytokeratin 5 Rabbit mAb                                                                                        |
| Remarks/Alias    | K5; KRT5; DDD; DDD1; EBS1; EBS2; EBS2A; EBS2B; EBS2C; EBS2D; EBS2E; EBS2F; KRT5A; Cytokeratin 5 (CK5)           |
| Species          | Human                                                                                                           |
| GeneID (Human)   | 3852                                                                                                            |
| GeneID           | 3852                                                                                                            |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 491-590 of human Cytokeratin 5 (KRT5) (P13647) |
| Source           | Rabbit                                                                                                          |
| Category         | Monoclonal Antibodies                                                                                           |
| Application      | WB, IF-P, IHC-P, ELISA                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | P13647                                                                                                          |
| Protein Weight   | 62kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |