


Anti- δ -Catenin/p120 Catenin Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - O60716

Item Number: CM0002842



Introduction

Rabbit monoclonal antibody against δ -Catenin/p120 Catenin (O60716). Validated for WB, IHC-P, IP, ELISA in human, mouse, rat. Regulates cadherin stability, transcription, and Rho GTPase signaling. Ideal for cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	δ -Catenin/p120 Catenin Rabbit mAb
Remarks/Alias	CAS, p120, BCD52, CTNND, P120CAS, P120CTN, p120(CAS), p120(CTN), δ -Catenin/p120 Catenin
Species	Human
GeneID (Human)	1500
GeneID	1500
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 869-968 of human δ -Catenin/p120 Catenin (O60716)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60716
Protein Weight	108kDa
Shipping	Ice bag