


Anti- Δ -Catenin/P120 Catenin Rabbit Polyclonal Antibody For Wb And Elisa - O60716

Item Number: CM0024613



Introduction

Rabbit polyclonal antibody against human δ -Catenin/p120 Catenin, validated for WB and ELISA. Reactive with human, mouse, rat. For cell adhesion and cancer research. Predicted MW 108 kDa.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] δ -Catenin/p120 Catenin Rabbit pAb
Remarks/Alias	CAS; p120; BCDS2; CTNND; P120CAS; P120CTN; p120(CAS); p120(CTN); in
Species	Human
Gene ID (Human)	1500
Gene ID	1500
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60716
Protein Weight	108kDa
Shipping	Ice bag