

Anti-Dap Kinase 1 (Dapk1) Polyclonal Antibody For Wb, Ihc-P, Elisa - P53355

Item Number: CM0027435



Introduction

High-quality rabbit polyclonal antibody targeting human DAP Kinase 1 (DAPK1), validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for apoptosis and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	DAP Kinase 1 (DAPK1) Rabbit pAb
Remarks/Alias	DAPK; ROCO3; DAP Kinase 1 (DAPK1)
Species	Human
Gene ID (Human)	1612
Gene ID	1612
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1141-1430 of human DAP Kinase 1 (DAPK1) (NP_004929.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53355
Protein Weight	160kDa
Shipping	Ice bag