

Anti-Dbf4 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Ubu7

Item Number: CM0002257



Introduction

High-quality anti-DBF4 rabbit monoclonal antibody for Western blot and ELISA. Targets human DBF4 (Q9UBU7), a key regulatory subunit of the CDC7 kinase involved in DNA replication and S phase progression. Suitable for cancer research.

[Learn More](#)

Parameter	Value
Product Name	DBF4 Rabbit mAb
Remarks/Alias	ASK; CHIF; DBF4A; ZDBF1; DBF4
Species	Human
Gene ID (Human)	10926
Gene ID	10926
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DBF4 (Q9UBU7).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UBU7
Protein Weight	77kDa
Shipping	Ice bag