


Anti-Dcaf13 Polyclonal Antibody For Wb, Elisa - Q9Nv06

Item Number: CM0013503



Introduction

Rabbit polyclonal antibody against human DCAF13, validated for WB and ELISA. Cross-reacts with mouse, rat. Recombinant immunogen (aa 298-597). Suitable for ribosome biogenesis, ubiquitination, and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	DCAF13 Rabbit pAb
Remarks/Alias	GM83; Sof1; WDSOF1; HSPC064; DCAF13
Species	Human
GeneID (Human)	25879
GeneID	25879
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 298-597 of human DCAF13 (NP_056235.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NV06
Protein Weight	51kDa
Shipping	Ice bag