

Anti-Dcaf17 Rabbit Polyclonal Antibody For Wb, Elisa - Q5H9S7

Item Number: CM0022974

Introduction

Rabbit polyclonal antibody against DCAF17 (C2orf37, C20orf37), a substrate receptor for CUL4-DDB1 E3 ubiquitin ligase. Validated for WB and ELISA, reacts with human and rat. Ideal for ubiquitination and disease research.

[Learn More](#)

Parameter	Value
Product Name	DCAF17 Rabbit pAb
Remarks/Alias	C2orf37; C20orf37; DCAF17
Species	Human
Gene ID (Human)	80067
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 235-335 of human DCAF17 (NP_079276.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q5H9S7
Protein Weight	59kDa
Shipping	Ice bag