

Anti-Dcdc2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uhg0

Item Number: CM0019674



Introduction

DCDC2 Rabbit pAb validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Designed for Wnt signaling, neuronal migration, and ciliogenesis research in neuroscience and cell biology.

[Learn More](#)

Parameter	Value
Product Name	DCDC2 Rabbit pAb
Remarks/Alias	NSC; RU2; RU2S; DCDC2A; DFNB66; NPHP19; DCDC2
Species	Human
Gene ID	51473
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human DCDC2 (NP_001182539)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHG0
Protein Weight	53kDa
Shipping	Ice bag