


Anti-Dcp1A Rabbit Polyclonal Antibody For Wb, Elisa - Q9Npi6

Item Number: CM0029338



Introduction

Rabbit polyclonal antibody against human DCP1A (SMIF), suitable for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for mRNA decay and TGF-beta signaling research.

[Learn More](#)

Parameter	Value
Product Name	DCP1A Rabbit pAb
Remarks/Alias	SMIF; SMAD4IP1; HSA275986; Nbla00360; DCP1A
Species	Human
GeneID (Human)	55802
GeneID	55802
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 303-582 of human DCP1A (NP_060873.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPI6
Protein Weight	63kDa
Shipping	Ice bag