

Anti-Dctd Polyclonal Antibody For Wb, Ihc-P, Elisa - P32321

Item Number: CM0028486



Introduction

Rabbit polyclonal antibody against human DCTD (dCMP deaminase), suitable for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, rat. Targets key enzyme in nucleotide metabolism and epigenetic salvage.

[Learn More](#)

Parameter	Specification
Product Name	DCTD Rabbit pAb
Remarks/Alias	DCTD
Species	Human
GeneID (Human)	1635
GeneID	1635
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of human DCTD (NP_001912.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32321
Protein Weight	20kDa
Shipping	Ice bag