


Anti-Ddah1 Rabbit Monoclonal Antibody For Wb And Elisa - O94760

Item Number: CM0004242



Introduction

High-quality rabbit monoclonal antibody targeting human DDAH1, validated in Western blotting and ELISA. Detects DDAH1 in human and rat samples. Essential for nitric oxide regulation and dimethylarginine metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	DDAH1 Rabbit mAb
Remarks/Alias	DDAH; DDAHI; DDAH-1; HEL-S-16; DDAH1
Species	Human
Gene ID (Human)	23576
Gene ID	23576
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 186-285 of human DDAH1 (O94760).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O94760
Protein Weight	31kDa
Shipping	Ice bag