

Anti-Ddah2 Rabbit Monoclonal Antibody For Wb, Elisa - 095865

Item Number: CM0012297

Introduction

DDAH2 rabbit monoclonal antibody validated for Western blot (WB) and ELISA applications. Detects human, mouse, and rat DDAH2, a protein involved in VEGF induction, insulin secretion, and mitochondrial dynamics. Ideal for vascular and immunology research.

[Learn More](#)

Parameter	Value
Product Name	DDAH2 Rabbit mAb
Remarks/Alias	G6a; DDAH; NG30; DDAHII; HEL-S-277; DDAH2
Species	Human
Gene ID (Human)	23564
Gene ID	23564
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DDAH2 (095865).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	095865
Protein Weight	30kDa
Shipping	Ice bag