

## Anti-Ddb1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q16531

Item Number: CM0025059



### Introduction

DDB1 rabbit polyclonal antibody for DNA damage repair and ubiquitination research. Validated in WB, IHC-P, IP, ELISA with human, mouse, rat reactivity. Ideal for studying UV-DDB complex and DCX E3 ligases.

[Learn More](#)

| Parameter        | Value                                                                                                     |
|------------------|-----------------------------------------------------------------------------------------------------------|
| Product Name     | DDB1 Rabbit pAb                                                                                           |
| Remarks/Alias    | XPE; DDBA; XAP1; XPCE; XPE-BF; UV-DDB1; WHIKERS; DDB1                                                     |
| Species          | Human                                                                                                     |
| GeneID (Human)   | 1642                                                                                                      |
| GeneID           | 1642                                                                                                      |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 1000-1140 of human DDB1 (NP_001914.3). |
| Source           | Rabbit                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                     |
| Application      | WB, IHC-P, IP, ELISA                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                         |
| SWISS            | Q16531                                                                                                    |
| Protein Weight   | 127kDa                                                                                                    |
| Shipping         | Ice bag                                                                                                   |