

Anti-Ddx11 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Fc9

Item Number: CM0029861

Introduction

Rabbit polyclonal antibody targeting human ATP-dependent DNA helicase DDX11, validated for Western blot and ELISA. Suitable for DNA replication, repair, and chromatin biology research.

[Learn More](#)

Parameter	Value
Product Name	DDX11 Rabbit pAb
Remarks/Alias	CHL1; KRG2; WABS; CHLR1; DDX11
Species	Human
GeneID (human)	1663
GeneID	1663
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human DDX11 (NP_689651.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96FC9
Protein Weight	108kDa
Shipping	Ice bag