


Anti-Ddx20 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uhi6

Item Number: CM0029492



Introduction

Rabbit polyclonal antibody against human DDX20 (Gemin3), a DEAD-box RNA helicase involved in snRNP assembly. Suitable for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for splicing research.

[Learn More](#)

Parameter	Value
Product Name	DDX20 Rabbit pAb
Remarks/Alias	DP103; GEMIN3; DDX20
Species	Human
GeneID (human)	11218
GeneID	11218
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-280 of human DDX20 (NP_009135.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHI6
Protein Weight	92kDa
Shipping	Ice bag