

Anti-Decorin Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P07585

Item Number: CM0019820

Introduction

Rabbit polyclonal antibody against human Decorin, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Protein weight 40kDa. Suitable for extracellular matrix and proteoglycan research.

[Learn More](#)

Parameter	Value
Product Name	Decorin Rabbit pAb
Remarks/Alias	CSCD, PG40, PGII, PGS2, DSPG2, SLRR1B, Decorin
Species	Human
GeneID (Human)	1634
GeneID	1634
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-359 of Decorin (NP_001911.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07585
Protein Weight	40kDa
Shipping	Ice bag