


Anti-Defa1 Polyclonal Antibody For Wb, If/Icc, Elisa - P59665

Item Number: CM0027109



Introduction

Rabbit polyclonal antibody against human neutrophil defensin 1 (DEFA1/HNP-1). Suitable for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Designed for innate immunity and antimicrobial defense studies.

[Learn More](#)

Parameter	Value
Product Name	DEFA1 Rabbit pAb
Remarks/Alias	HP1; MRS; DEF1; HP-1; DEFA2; HNP-1; DEFA1
Species	Human
GeneID (Human)	1667
GeneID	1667
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-94 of human DEFA1 (NP_004075.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P59665
Protein Weight	10kDa
Shipping	Ice bag