

Anti-Desmin Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P17661

Item Number: CM0011567

Introduction

Anti-Desmin Rabbit Monoclonal Antibody validated for WB, IF-P, IHC-P, ELISA in human, mouse, rat. Desmin is a type III intermediate filament essential for muscle structure. Ideal for skeletal and cardiac muscle research.

[Learn More](#)

Parameter	Value
Product Name	Desmin Rabbit mAb
Remarks/Alias	CSM1; CSM2; CDCD3; LGMD1D; LGMD1E; LGMD2R; Desmin
Species	Human
GeneID (Human)	1674
GeneID	1674
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 371-470 of human Desmin (P17661).)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17661
Protein Weight	54 kDa
Shipping	Ice bag