

Anti-Desmoglein 3 Rabbit Monoclonal Antibody For Wb And Elisa - P32926

Item Number: CM0004091



Introduction

Desmoglein-3 rabbit monoclonal antibody for WB and ELISA. Recognizes human Desmoglein-3 (P32926), a desmosomal cadherin key for keratinocyte adhesion. Ideal for studies of pemphigus, cell junctions, and epithelial barrier research.

[Learn More](#)

Parameter	Value
Product Name	Desmoglein 3 Rabbit mAb
Remarks/Alias	PVA; ABOLM; CDHF6; Desmoglein 3
Species	Human
Gene ID (Human)	1830
Gene ID	1830
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 900-999 of human Desmoglein 3 (P32926).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P32926
Protein Weight	108 kDa
Shipping	Ice bag