

Anti-Det1 Polyclonal Antibody For Wb And Elisa - P48732

Item Number: CM0019278

Introduction

This rabbit polyclonal antibody targets Arabidopsis thaliana DET1, a key regulator of light signal transduction and photomorphogenesis. Validated for WB and ELISA, it's ideal for plant developmental studies and ubiquitination pathway research.

[Learn More](#)

Parameter	Value
Product Name	DET1 Rabbit pAb
Remarks/Alias	ATDET1; DE-ETIOLATED 1; FUS2; FUSCA 2; T9A4.17; DET1
Species	Arabidopsis thaliana
GeneID	826609
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 350-450 of Arabidopsis thaliana DET1 (NP_192756.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	P48732
Protein Weight	62kDa
Shipping	Ice bag