

Anti-Dgka Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P23743

Item Number: CM0016205



Introduction

Rabbit polyclonal antibody against DGKA, validated for WB, IHC-P, and ELISA. Recognizes human, mouse, and rat diacylglycerol kinase alpha, a key enzyme converting DAG to PA. Ideal for lipid signaling and immune cell research.

[Learn More](#)

Parameter	Value
Product Name	DGKA Rabbit pAb
Remarks/Alias	DAGK; DAGK1; DGK-alpha; DGKA
Species	Human
Gene ID (Human)	1606
Gene ID	1606
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human DGKA (NP_001336.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23743
Protein Weight	83 kDa
Shipping	Ice bag