


Anti-Dlat Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P10515

Item Number: CM0014455



Introduction

DLAT Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Targets the 69kDa dihydrolipoamide acetyltransferase (PDC-E2), a key enzyme in the pyruvate dehydrogenase complex and primary biliary cirrhosis autoantigen.

[Learn More](#)

Parameter	Value
Product Name	DLAT Rabbit pAb
Remarks/Alias	E2; PBC; DLTA; PDCE2; PDC-E2; DLAT
Species	Human
GeneID	1737
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 87-270 of human DLAT (NP_001922.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10515
Protein Weight	69kDa
Shipping	Ice bag